



# Anti-SLC35B4 polyclonal antibody (CABT-B8938)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A peptide from the intermediate residues of human SLC35B4 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Protein A affinity chromatography.
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Preservative	0.02% Sodium Azide
Storage	The antibody is stable for at least one year when stored at -20°C after reconstitution. Avoid repeated freeze/thaw cycles.

## BACKGROUND

Introduction	Glycosyltransferases, such as SLC35B4, transport nucleotide sugars from the cytoplasm where they are synthesized.
--------------	---

## GENE INFORMATION

Entrez Gene ID [84912](#)

---

UniProt ID [Q969S0](#)

---