



Anti-TREML4 Polyclonal Antibody (CABT-B8944)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein. Sequence: TCCFHLCCCCSWPQGAVPEELHKHPGQTL LLQCQYSPKRGPYQPKSWCQQTSPSRCTL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

GENE INFORMATION

Entrez Gene ID	224840
UniProt ID	Q3LRV9