



Anti-YTHDF1 polyclonal antibody (CABT-B8963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	YTHDF1 Antibody detects endogenous levels of YTHDF1
Immunogen	A synthesized peptide derived from human YTHDF1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C.Stable for 12 months from date of receipt