



Anti-PRRC2B polyclonal antibody (CABT-B8985)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human PRRC2B
Immunogen	Recombinant protein.Sequence: EREESTLKKGDCRDSWRSNKGCSGDHSGLDKSRGPRAFGRALPPRLSNCGYGRRTFVSK ESPHWQSKSPGSSWQEYGPSDTCGSRRPTRDRDYVPDSYRHPDAFGGR
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywordsproline rich coiled-coil 2B; BAT2L; BAT2L1; LQFBS-1; KIAA0515; PRRC2B

GENE INFORMATION

Entrez Gene ID[84726](#)

UniProt ID[Q5JSZ5](#)
