



Anti-PRRC2C polyclonal antibody (CABT-B8986)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human PRRC2C
Immunogen	Recombinant protein. Sequence: DVRPHHTDANNQSACFEAPDQKTLSTPQEERISAVESQPSRKRSVSHGSNHTQKPDEQRS EPSAGIPKVTSRCIDSKEP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords

proline rich coiled-coil 2C; PRRC2C; XTP2; BAT2D1; BAT2L2; BAT2-iso

GENE INFORMATION

Entrez Gene ID[23215](#)**UniProt ID**[Q9Y520](#)
