



Rabbit Anti-T7 RNA Polymerase polyclonal antibody [Biotin] (CABT-B8990B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Rabbit polyclonal Antibody specifically binds to T7 RNA Polymerase.
Immunogen	Recombinant T7 RNA Polymerase (full-length)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Affinity purified.
Conjugate	Biotin
Applications	ELISA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	100 µl, 500 µl
Buffer	PBS with 0.2M Glycine
Preservative	None
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Introduction

T7 RNA Polymerase is an RNA polymerase from the T7 bacteriophage that catalyzes the formation of RNA in the 5'→ 3' direction.

Keywords

T7 RNA Polymerase; RNA Polymerase
