



# Rabbit Anti-T7 RNA Polymerase polyclonal antibody[HRP] (CABT-B8990H)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit polyclonal Antibody specifically binds to T7 RNA Polymerase, HRP conjugated
<b>Immunogen</b>	Recombinant T7 RNA Polymerase (full-length)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Purification</b>	The pAb was purified by protein A affinity chromatography, ultrafiltration and then subjected to 0.2 micron sterile filtration to get the bulk of high purity
<b>Conjugate</b>	HRP
<b>Applications</b>	ELISA, IF
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µg
<b>Buffer</b>	PBS 7.4 with 30% glycerol, Casein/BSA, and Proclin 300
<b>Storage</b>	Store -20°C for long term.
<b>Ship</b>	Wet ice

## BACKGROUND

**Introduction**

T7 RNA Polymerase is an RNA polymerase from the T7 bacteriophage that catalyzes the formation of RNA in the 5'→ 3' direction

---

**Keywords**

T7 RNA Polymerase

---