



Rabbit Anti-Human OR4F6 Polyclonal Antibody (CABT-B9133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to Olfactory receptor 4F6.
Specificity	This antibody detects endogenous levels of Olfactory receptor 4F6 protein.
Target	Olfactory receptor 4F6
Immunogen	Synthesized peptide derived from human Olfactory receptor 4F6 (Internal)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, IF, ELISA
Molecular Weight	35 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS containing 50% glycerol and 0.5% BSA

Preservative	0.02% Sodium Azide
---------------------	--------------------

Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.
----------------	--

BACKGROUND

Introduction	OR4F6 belongs to the G-protein coupled receptor 1 family and is an odorant receptor.
---------------------	--

Keywords	OR4F6;olfactory receptor, family 4, subfamily F, member 6;olfactory receptor, family 4, subfamily F, member 12;OR4F12;OR4F12;
-----------------	---

GENE INFORMATION

Entrez Gene ID	390648
-----------------------	------------------------

UniProt ID	Q8NGB9
-------------------	------------------------