



Rabbit Anti-Human ZNF682 Polyclonal Antibody (CABT-B9150)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to ZNF682.
Specificity	This antibody detects endogenous levels of ZNF682 protein.
Target	ZNF682
Immunogen	Synthesized peptide derived from human ZNF682 (N-terminal)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	58 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl
Buffer	PBS containing 50% glycerol and 0.5% BSA

Preservative	0.02% Sodium Azide
---------------------	--------------------

Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.
----------------	--

BACKGROUND

Introduction	May be involved in transcriptional regulation.
---------------------	--

Keywords	ZNF682;zinc finger protein 682;BC39498_3;
-----------------	---

GENE INFORMATION

Entrez Gene ID	91120
-----------------------	-----------------------

UniProt ID	O95780
-------------------	------------------------