



Rabbit Anti-Human CSPP1 Polyclonal Antibody (CABT-B9156)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to CSPP1
Target	CSPP1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	138 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS containing 50% glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction

CSPP1 interacts with centrosomes and microtubules and may play a role in cell-cycle-dependent microtubule organization.

Keywords

CSPP1;centrosome and spindle pole associated protein 1;centrosome and spindle pole-associated protein 1;CSPP;FLJ22490;FLJ38886;

GENE INFORMATION

Entrez Gene ID

[79848](#)

UniProt ID

[Q1MSJ5](#)