



Rabbit Anti-Human OR1L1 Polyclonal Antibody (CABT-B9163)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to OR1L1
Target	OR1L1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	39 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS containing 50% glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction	Odorant receptor. Belongs to the G-protein coupled receptor 1 family.
Keywords	OR1L1;olfactory receptor, family 1, subfamily L, member 1;HG23;OR1L2;OR9-C;OR9-27;olfactory receptor 1L1;olfactory receptor 1L2;olfactory receptor 9-C;olfactory receptor OR9-27;olfactory receptor, family 1, subfamily L, member 2;

GENE INFORMATION

Entrez Gene ID	26737
UniProt ID	Q8NH94