



Rabbit Anti-Human OR1L4 Polyclonal Antibody (CABT-B9164)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to OR1L4
Target	OR1L4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	34 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS containing 50% glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction	Odorant receptor. Belongs to the G-protein coupled receptor 1 family.
Keywords	OR1L4;olfactory receptor, family 1, subfamily L, member 4;OR1L5;OR9-E;OR9-29;OST046;olfactory receptor 1L4;olfactory receptor 1L5;olfactory receptor 9-E;olfactory receptor OR9-29;olfactory receptor, family 1, subfamily L, member 5;

GENE INFORMATION

Entrez Gene ID	254973
UniProt ID	Q8NGR5