



# Rabbit Anti-Human OR2V1 Polyclonal Antibody (CABT-B9165)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Polyclonal to OR2V1
<b>Target</b>	OR2V1
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Molecular Weight</b>	34 kDa
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	100 µl, 200 µl
<b>Buffer</b>	PBS containing 50% glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C, and avoid repeat freeze-thaw cycles.

# BACKGROUND

<b>Introduction</b>	Odorant receptor. Belongs to the G-protein coupled receptor 1 family.
<b>Keywords</b>	OR2V1;olfactory receptor, family 2, subfamily V, member 1;OR2V1P;OST265;olfactory receptor 2V1;olfactory receptor, family 2, subfamily V, member 1 pseudogene;

# GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">26693</a>
<b>UniProt ID</b>	<a href="#">Q8NHB1</a>