



# Mouse anti-Human Maspin monoclonal antibody, clone H278-80 (CABT-B9233)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant Human Maspin
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	H278-80
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	0.5 mg/ml
Size	100 µg
Buffer	Aqueous buffered solution containing ≤0.09% sodium azide.
Storage	Store undiluted at 4°C.

## BACKGROUND

**Keywords** ITGB1BP2; integrin beta 1 binding protein (melusin) 2; integrin beta-1-binding protein 2;

