



Mouse anti-Rat Synapsin IIa monoclonal antibody, clone 2 (CABT-B9330)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Rat Synapsin IIa aa. 451-586
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Rat, Human
Clone	2
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	Unconjugated
Applications	WB; IP; IHC; IF
Format	Liquid
Concentration	250 µg/ml
Size	50 µg
Buffer	Aqueous buffered solution containing BSA, glycerol, and ≤0.09% sodium azide.
Storage	Store undiluted at -20°C.

BACKGROUND

Keywords	SYN2; synapsin II; SYNII; synapsin-2;
-----------------	---------------------------------------

GENE INFORMATION

Entrez Gene ID [6854](#)

UniProt ID [A0A087X2E3](#)
