



Mouse anti-Human Tom20 monoclonal antibody, clone 30/Upn30 (CABT-B9342)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human Tom20 aa. 47-145
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human, Rat, Dog
Clone	30/Upn30
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	Unconjugated
Applications	WB; IF
Format	Liquid
Concentration	250 µg/ml
Size	50 µg
Buffer	Aqueous buffered solution containing BSA, glycerol, and ≤0.09% sodium azide.
Storage	Store undiluted at -20°C.

BACKGROUND

Keywords TOMM20; translocase of outer mitochondrial membrane 20 homolog (yeast); MAS20; MOM19;

TOM20; mitochondrial import receptor subunit TOM20 homolog; mitochondrial 20 kDa outer membrane protein; outer mitochondrial membrane receptor Tom20; translocase of outer mitochondrial membrane 20 homolog type II;

GENE INFORMATION

Entrez Gene ID	9804
----------------	----------------------

UniProt ID	A0A024R3W2
------------	----------------------------
