



Mouse anti-Rat Tomosyn monoclonal antibody, clone 26/Upnptzo (CABT-B9343)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Rat Tomosyn aa. 910-1105
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Rat, Human, Mouse
Clone	26/Upnptzo
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	Unconjugated
Applications	WB; Bioimaging; IF
Format	Liquid
Concentration	250 µg/ml
Size	50 µg
Buffer	Aqueous buffered solution containing BSA, glycerol, and ≤0.09% sodium azide.
Storage	Store undiluted at -20°C.

BACKGROUND

Keywords STXBP5; syntaxin binding protein 5 (tomosyn); LGL3; AW742610; mKIAA1006;

0710001E20Rik; 4930565N16Rik; syntaxin-binding protein 5; tomosyn 1; tomosyn-1; lethal(2)
giant larvae protein homolog 3;

GENE INFORMATION

Entrez Gene ID [78808](#)

UniProt ID [D3Z079](#)
