



Mouse anti-Human Underphosphorylated Retinoblastoma Protein (RB) monoclonal antibody, clone H00-650 (CABT-B9356)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant full-length human Rb
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human, Mouse
Clone	H00-650
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	Unconjugated
Applications	IP
Format	Liquid
Concentration	0.5 mg/ml
Size	100 µg
Buffer	Aqueous buffered solution containing ≤0.09% sodium azide.
Storage	Store undiluted at 4°C.

BACKGROUND

Keywords

RB1; retinoblastoma 1; RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; retinoblastoma-associated protein; retinoblastoma susceptibility protein; prepro-retinoblastoma-associated protein; protein phosphatase 1, regulatory subunit 130;

GENE INFORMATION

Entrez Gene ID[5925](#)**UniProt ID**[A0A024RDV3](#)