



# Mouse anti-Human WIP monoclonal antibody, clone L43-935 [PE] (CABT-B9360)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Phosphorylated Human WIP Peptide
Isotype	IgG1, $\kappa$
Source/Host	Mouse
Species Reactivity	Human
Clone	L43-935
Purification	The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.
Conjugate	PE
Applications	FC
Format	Liquid
Size	50 tests
Buffer	Aqueous buffered solution containing BSA and $\leq 0.09\%$ sodium azide.
Storage	Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

## BACKGROUND

Keywords	WIPF1; WAS/WASL interacting protein family, member 1; WIP; PRPL-2; WASPIP; WAS/WASL-interacting protein family member 1; protein PRPL-2; WASP-interacting protein; Wiskott-Aldrich syndrome protein interacting protein;
----------	--

# GENE INFORMATION

Entrez Gene ID [7456](#)

---

UniProt ID [O43516](#)

---