



Mouse anti-Human ADCK1 monoclonal antibody, clone 2F9 (CABT-B9375)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ADCK1 (AAH58906.1, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2F9
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords ADCK1; aarF domain containing kinase 1; uncharacterized aarF domain-containing protein kinase 1;

GENE INFORMATION

Entrez Gene ID [57143](#)

UniProt ID [Q86TW2](#)
