



Mouse anti-Human AHSA2 monoclonal antibody, clone 4H22 (CABT-B9376)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	AHSA2 (NP_689605.1, 1 a.a. ~ 137 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4H22
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords AHSA2; AHA1, activator of heat shock 90kDa protein ATPase homolog 2 (yeast); Hch1;

activator of 90 kDa heat shock protein ATPase homolog 2;

GENE INFORMATION

Entrez Gene ID [130872](#)

UniProt ID [Q719I0](#)
