



Rat anti-Mouse AMIGO3 monoclonal antibody, clone 482402 (CABT-B9377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse AMIGO3 Thr20-Pro378
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Mouse
Clone	482402
Purification	Protein A/G
Conjugate	Unconjugated
Applications	IHC-Fr, WB
Format	Lyophilized
Concentration	0.5 mg/ml
Size	100 µg
Buffer	PBS with 5% trehalose
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

AMIGO3; adhesion molecule with Ig like domain 3; ali3; mKIAA1851; E430002N15Rik; amphoterin-induced protein 3; amigo3; AMIGO-3; alivin 3; alivin-3; amphoterin induced gene and ORF 3;

GENE INFORMATION

Entrez Gene ID[320844](#)**UniProt ID**[Q8C2S7](#)