



# Mouse anti-Human ARID5A monoclonal antibody, clone Q032 (CABT-B9379)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Recombinant fragment corresponding to a region within amino acids 1 and 259 of ARID5A
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	Q032
Purification	Protein G
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	1 mg/ml
Size	100 µl
Buffer	PBS
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction** Members of the ARID protein family, including ARID5A, have diverse functions but all appear to play important roles in development, tissue-specific gene expression, and regulation of cell growth (Patsialou et al., 2005 [PubMed 15640446]).[supplied by OMIM, Mar 2008]

**Keywords** ARID5A; AT rich interactive domain 5A (MRF1-like); MRF1; MRF-1; RP11-363D14; AT-rich interactive domain-containing protein 5A; RFVG5814; modulator recognition factor 1; modulator recognition factor I; ARID domain-containing protein 5A;

## GENE INFORMATION

**Entrez Gene ID** [10865](#)

**UniProt ID** [Q03989](#)