



Mouse anti-Human ARM6 monoclonal antibody, clone 3I9 (CABT-B9380)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ARM6 (AAH67125, 1 a.a. ~ 476 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3I9
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, IM, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

The function of this genes protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010]

Keywords

ARMC6; armadillo repeat containing 6; R30923_1; armadillo repeat-containing protein 6;

GENE INFORMATION

Entrez Gene ID

[93436](#)

UniProt ID

[Q6NXE6](#)