



Mouse anti-Human BLZF1 monoclonal antibody, clone 3F20 (CABT-B9384)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	BLZF1 (NP_003657, 301 a.a. ~ 400 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3F20
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords BLZF1; basic leucine zipper nuclear factor 1; JEM1; JEM-1; JEM-1s; GOLGIN-45; Golgin-45;
JEM-1short protein; cytoplasmic protein; p45 basic leucine-zipper nuclear factor;

GENE INFORMATION

Entrez Gene ID [8548](#)

UniProt ID [A0A024R921](#)
