



Mouse anti-Human C10orf54 monoclonal antibody, clone DM4086 (CABT-B9387)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Human C10orf54
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	DM4086
Purification	Protein A
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2, with 40% glycerol
Preservative	0.02% sodium azide
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Keywords	C10ORF54; chromosome 10 open reading frame 54; B7H5; GI24; B7-H5; SISP1; PP2135; platelet receptor Gi24; sisp-1; stress induced secreted protein 1; stress-induced secreted
----------	---

GENE INFORMATION

Entrez Gene ID	64115
----------------	-----------------------

UniProt ID	Q9H7M9
------------	------------------------
