



Mouse anti-Human C1r monoclonal antibody, clone 370205 (CABT-B9388)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Mouse myeloma cell line NS0-derived recombinant human Complement Component C1r Ser18-Asp705
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	370205
Purification	Protein A/G
Conjugate	Unconjugated
Applications	IM, WB
Format	Lyophilized
Concentration	0.5 mg/ml
Size	100 µg
Buffer	PBS with 5% trehalose
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

C1R; complement component 1, r subcomponent; complement C1r subcomponent;

GENE INFORMATION

Entrez Gene ID[715](#)**UniProt ID**[P00736](#)
