



Mouse anti-Human CBX8 monoclonal antibody [Biotin] (CABT-B9390)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	1-70 of 389 of CBX8 Protein.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Purification	purified
Conjugate	Biotin
Applications	IM
Format	Liquid
Concentration	0.5 mg/ml
Size	100 µg
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Maintain refrigerated at 2-8°C for up to 1 month. For long term storage store at -20°C

BACKGROUND

Keywords CBX8; chromobox homolog 8; PC3; RC1; chromobox protein homolog 8; polycomb 3; rectachrome 1; Pc class 3 homolog; chromobox homolog 8 (Pc class homolog, Drosophila);

GENE INFORMATION

Entrez Gene ID [57332](#)

UniProt ID [Q9HC52](#)
