



Mouse anti-Human CCDC83 monoclonal antibody, clone 4F4 (CABT-B9391)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CCDC83 (NP_775827, 341 a.a. ~ 444 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4F4
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords CCDC83; coiled-coil domain containing 83; HSD9; CT148; coiled-coil domain-containing protein 83;

GENE INFORMATION

Entrez Gene ID [220047](#)

UniProt ID [Q8IWF9](#)
