



Mouse anti-Human CENPV monoclonal antibody, clone 2E0 (CABT-B9393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CENPV (NP_859067.2, 173 a.a. ~ 272 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2E0
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords CENPV; centromere protein V; p30; PRR6; CENP-V; 3110013H01Rik; proline rich 6; nuclear

GENE INFORMATION

Entrez Gene ID [201161](#)

UniProt ID [Q7Z7K6](#)
