



Mouse anti-Human COMMD3 monoclonal antibody, clone 3F3 (CABT-B9400)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	COMMD3 (AAH22898, 1 a.a. ~ 195 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3F3
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords COMMD3; COMM domain containing 3; BUP; C10orf8; COMM domain-containing protein 3; protein Bup;

GENE INFORMATION

Entrez Gene ID [23412](#)

UniProt ID [Q9UBI1](#)
