



Mouse anti-Human DEPDC7 monoclonal antibody, clone 5E0 (CABT-B9405)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DEPDC7 (NP_631899, 393 a.a. ~ 499 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5E0
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, IM, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords DEPDC7; DEP domain containing 7; TR2; dJ85M6.4; DEP domain-containing protein 7; novel 58.3 KDA protein; dJ85M6.4 (novel 58.3 KDA protein);

GENE INFORMATION

Entrez Gene ID [91614](#)

UniProt ID [Q96QD5](#)
