



Mouse anti-Human DMRTB1 monoclonal antibody, clone 6F8 (CABT-B9408)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DMRTB1 (AAH29566.1, 1 a.a. ~ 189 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	6F8
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords DMRTB1; DMRT-like family B with proline-rich C-terminal, 1; doublesex- and mab-3-related

transcription factor B1;

GENE INFORMATION

Entrez Gene ID [63948](#)

UniProt ID [Q96MA1](#)
