



Mouse anti-Human DNAJC19 monoclonal antibody, clone 4I5 (CABT-B9409)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DNAJC19 (NP_660304, 20 a.a. ~ 116 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4I5
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

DNAJC19; DnaJ (Hsp40) homolog, subfamily C, member 19; Tim14; AA959924;
1810055D05Rik; mitochondrial import inner membrane translocase subunit TIM14;

GENE INFORMATION

Entrez Gene ID[67713](#)**UniProt ID**[G3X8S1](#)