



Mouse anti-Human DNAJC27 monoclonal antibody, clone 4D2 (CABT-B9410)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DNAJC27 (AAH34049, 1 a.a. ~ 273 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4D2
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

DNAJC27; DnaJ (Hsp40) homolog, subfamily C, member 27; RBJ; RabJS; dnaJ homolog subfamily C member 27; Ras-associated protein Rap1; rab and DnaJ domain-containing protein;

GENE INFORMATION

Entrez Gene ID[51277](#)**UniProt ID**[Q9NZQ0](#)
