



Mouse anti-Human DPCR1 monoclonal antibody, clone 324302 (CABT-B9411)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | NS0 mouse myeloma cell line transfected with human DPCR-1 Met1-Arg235 |
| Isotype | IgG2a |
| Source/Host | Mouse |
| Species Reactivity | Human |
| Clone | 324302 |
| Purification | Protein A/G |
| Conjugate | Unconjugated |
| Applications | FC |
| Format | Lyophilized |
| Concentration | 0.5 mg/ml |
| Size | 100 µg |
| Buffer | PBS with 5% trehalose |
| Preservative | no preservative |
| Storage | -20°C, Avoid Freeze/Thaw Cycles |

BACKGROUND

Keywords

DPCR1; diffuse panbronchiolitis critical region 1; PBLT; diffuse panbronchiolitis critical region protein 1;

GENE INFORMATION

Entrez Gene ID

[135656](#)

UniProt ID

[E9PEI6](#)
