



Mouse anti-Human EFHC1 monoclonal antibody, clone 5F8 (CABT-B9416)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	EFHC1 (NP_060570.1, 270 a.a. ~ 378 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5F8
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords EFHC1; EF-hand domain (C-terminal) containing 1; mRib72-1; myoclonin1; 1700029F22Rik;

GENE INFORMATION

Entrez Gene ID [71877](#)

UniProt ID [B2CKC6](#)
