



Mouse anti-Human eIF3i monoclonal antibody, clone 4C20 (CABT-B9418)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	EIF3I (NP_003748.1, 226 a.a. ~ 325 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4C20
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

EIF3I; eukaryotic translation initiation factor 3, subunit I; TRIP1; EIF3S2; TRIP-1; PRO2242; eIF3-p36; eIF3-beta; eukaryotic translation initiation factor 3 subunit I; eIF3 p36; eIF-3-beta; predicted protein of HQ2242; TGFbeta receptor-interacting protein 1; TGF-beta receptor-interacting protein 1; eukaryotic translation initiation factor 3 subunit 2; eukaryotic translation initiation factor 3, subunit 2 beta, 36kDa; eukaryotic translation initiation factor 3, subunit 2 (beta, 36kD);

GENE INFORMATION

Entrez Gene ID[8668](#)**UniProt ID**[Q13347](#)