



# Mouse anti-Human ELTD1 monoclonal antibody, clone DM5275 (CABT-B9419)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Human ELTD1
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	DM5275
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	PBS, pH 7.2, with 40% glycerol
<b>Preservative</b>	0.02% sodium azide
<b>Storage</b>	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

## BACKGROUND

<b>Keywords</b>	ELTD1; EGF, latrophilin and seven transmembrane domain containing 1; ETL; KPG_003; EGF, latrophilin and seven transmembrane domain-containing protein 1; EGF-TM7-latrophilin-related
-----------------	--

protein;

---

## GENE INFORMATION

---

Entrez Gene ID	<a href="#">64123</a>
----------------	-----------------------

---

UniProt ID	<a href="#">Q9HBW9</a>
------------	------------------------

---