



# Mouse anti-Human FAM174A monoclonal antibody, clone 2G4 (CABT-B9423)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	FAM174A (NP_940909.1, 1 a.a. ~ 190 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Isotype</b>	IgG2a, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	2G4
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	no preservative
<b>Storage</b>	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Keywords** FAM174A; family with sequence similarity 174, member A; TMEM157; UNQ1912;

HGS\_RE408; membrane protein FAM174A; transmembrane protein 157; HCV NS5A-transactivated protein 6; hepatitis C virus NS5A-transactivated protein 6;

---

## GENE INFORMATION

Entrez Gene ID [345757](#)

---

UniProt ID [Q8TBP5](#)

---