



Mouse anti-Human FBXL22 monoclonal antibody, clone 2C7 (CABT-B9427)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FBXL22 (NP_976307, 1 a.a. ~ 110 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2C7
Purification	Protein A
Conjugate	Unconjugated
Applications	ICC, IF
Format	Liquid
Size	50 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction This gene encodes a member of the F-box protein family. This F-box protein interacts with S-

phase kinase-associated protein 1A and cullin in order to form SCF complexes which function as ubiquitin ligases.[provided by RefSeq, Sep 2010]

Keywords	FBXL22; F-box and leucine-rich repeat protein 22; Fbl22; F-box and leucine-rich protein 22; F-box/LRR-repeat protein 22;
-----------------	--

GENE INFORMATION

Entrez Gene ID	283807
-----------------------	------------------------

UniProt ID	Q6P050
-------------------	------------------------
