



Mouse anti-Human GATA5 monoclonal antibody, clone 2C22 (CABT-B9441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Purified recombinant fragment of human GATA5 expressed in E. Coli.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2C22
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	ascites
Preservative	0.03% sodium azide
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Keywords GATA5; GATA binding protein 5; gta5; cb634; fau/gata5; transcription factor GATA-5; fau; faust; GATA-binding protein 5;

GENE INFORMATION

Entrez Gene ID [30482](#)

UniProt ID [Q9W6U0](#)
