



Mouse anti-Human GDE1 monoclonal antibody, clone 3I7 (CABT-B9442)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GDE1 (NP_057725, 31 a.a. ~ 138 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3I7
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

GDE1; glycerophosphodiester phosphodiesterase 1; MIR16; 363E6.2; RGS16-interacting membrane protein; membrane interacting protein of RGS16; membrane-interacting protein of RGS16;

GENE INFORMATION

Entrez Gene ID[51573](#)**UniProt ID**[Q9NZC3](#)
