



Mouse anti-Human GIMAP2 monoclonal antibody, clone 2F20 (CABT-B9444)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GIMAP2 (NP_056475.1, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2F20
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008]

Keywords

GIMAP2; GTPase, IMAP family member 2; IAN12; IMAP2; HIMAP2; GTPase IMAP family member 2; immunity associated protein 2; immunity-associated protein 2; immune-associated nucleotide-binding protein 12;

GENE INFORMATION

Entrez Gene ID

[26157](#)

UniProt ID

[A0A090N8H4](#)