



Mouse anti-Human GPX7 monoclonal antibody, clone 3805 (CABT-B9448)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human NPGPx
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	3805
Purification	Protein G
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Concentration	1 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords GPX7; glutathione peroxidase 7; GPX6; CL683; GPx-7; NPGPx; GSHPx-7; glutathione peroxidase 6; non-selenocysteine containing phospholipid hydroperoxide glutathione peroxidase;

GENE INFORMATION

Entrez Gene ID [2882](#)

UniProt ID [Q96SL4](#)
