



Mouse anti-Human HCP5 monoclonal antibody, clone 4H8 (CABT-B9451)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HCP5 (NP_006665, 3 a.a. ~ 103 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4H8
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 μ g
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

HCP5; HLA complex P5 (non-protein coding); P5-1; 6S2650E; D6S2650E;

GENE INFORMATION

Entrez Gene ID

[10866](#)
