



# Mouse anti-Human KLHL22 monoclonal antibody, clone 8C20 (CABT-B9464)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Purified recombinant fragment of human KLHL22 expressed in E. Coli.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human, Mouse
Clone	8C20
Conjugate	Unconjugated
Applications	ICC, IF, WB
Format	Liquid
Size	100 µl
Buffer	ascites
Preservative	0.03% sodium azide
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

## BACKGROUND

**Keywords** KLHL22; kelch-like family member 22; KELCHL; kelch-like protein 22; kelch-like 22;

## GENE INFORMATION

Entrez Gene ID [84861](#)

---

UniProt ID [Q53GT1](#)

---