



# Mouse anti-Human LIN7A monoclonal antibody, clone Q980 (CABT-B9466)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant fragment corresponding to a region within amino acids 12 and 227 of LIN7A
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	Q980
<b>Purification</b>	Protein G
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC, IF, WB
<b>Format</b>	Liquid
<b>Concentration</b>	1 mg/ml
<b>Size</b>	100 µl
<b>Buffer</b>	PBS
<b>Preservative</b>	no preservative
<b>Storage</b>	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Keywords**

LIN7A; lin-7 homolog A (C. elegans); LIN7; VELI1; LIN-7A; MALS-1; TIP-33; protein lin-7 homolog A; hLin-7; veli-1; mammalian LIN-7 1; vertebrate LIN7 homolog 1; tax interaction protein 33; vertebrate lin-7 homolog 1; mammalian lin-seven protein 1;

---

## GENE INFORMATION

**Entrez Gene ID**[8825](#)**UniProt ID**[O14910](#)