



# Mouse anti-Human METTL2B monoclonal antibody, clone 3B0 (CABT-B9472)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	METTL2B (NP_060866.1, 41 a.a. ~ 108 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Isotype</b>	IgG2a, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	3B0
<b>Purification</b>	Affinity chromatography
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Concentration</b>	0.2-1.0 mg/ml
<b>Size</b>	100 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	no preservative
<b>Storage</b>	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

This gene is a member of a family of methyltransferases that share homology with, but are distinct from, the UbiE family of methyltransferases. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized. [provided by RefSeq, Jul 2008]

**Keywords**

METTL2B; methyltransferase like 2B; METL; METTL2; METTL2A; PSENIP1; methyltransferase-like protein 2B;

## GENE INFORMATION

**Entrez Gene ID**

[55798](#)

**UniProt ID**

[Q6P1Q9](#)