



Mouse anti-Human MPPE1 monoclonal antibody, clone 4F5 (CABT-B9475)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MPPE1 (AAH02877.1, 34 a.a. ~ 333 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4F5
Purification	Protein A
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords MPPE1; metallophosphoesterase 1; PGAP5; metallo phosphoesterase; post-GPI attachment to

GENE INFORMATION

Entrez Gene ID [65258](#)

UniProt ID [Q53F39](#)
